

# ATLANTIS™ abutments

## — Screw torque recommendations

Manufacturer	Implant	Torque Ncm
<b>BioHorizons</b>	Internal 3.5, 4.5, 5.7	30
	Tapered Internal 3.5, 4.5, 5.7	30
	External implant 3.5	30
	External implant 4.0, 5.0	30
<b>Biomet 3i</b>	Certain™ 3.25, 4.0, 4/3, 5.0, 5/4, 6.0	20
	Certain™ Prevail™ 3/4/3 (3.4), 4/5/4, 5/6/5, 4/3, 5/4, 6/5	20
	Certain™ XP 4/5, 5/6	20
	MicroMiniplant™ 3.25	20
	Miniplant™ 3.25	20
	Osseotite™ 3.75, 4.0, 5.0, 6.0	20
	Osseotite™ XP 4/5, 5/6	20
<b>Camlog</b>	Screw-line Implant 3.3	20
	Screw-line / Root-line Implant 3.8, 4.3, 5.0, 6.0	20
<b>DENTSPLY Implants</b>  Note: OsseoSpeed™ implants are part of the ASTRA TECH Implant System™	ANKYLOS® C/X 3.5, 4.5, 5.5, 7.0	15
	OsseoSpeed™ EV 3.0 – Green	25
	OsseoSpeed™ EV 3.6 – Purple	25
	OsseoSpeed™ EV 4.2 – Yellow	25
	OsseoSpeed™ EV 4.8 – Blue	25
	OsseoSpeed™ EV 5.4 – Brown	25
	OsseoSpeed™ / OsseoSpeed™ TX 3.0 – Yellow	15
	OsseoSpeed™ / OsseoSpeed™ TX 3.5, 4.0 – Aqua	20
	OsseoSpeed™ / OsseoSpeed™ TX 4.5, 5.0 – Lilac	25
	OsseoSpeed™ 4.0 - Blue (Previous)	25
	OsseoSpeed™ TX Profile 4.5, 5.0	25
	XIVE® S 3.0	24
	XIVE® S 3.4, 3.8, 4.5, 5.5	24
<b>Keystone Dental</b>	PrimaConnex™ SD 3.3, 3.5, Genesis™ 3.8	30
	PrimaConnex™ RD 4.0, 4.1, Genesis™ 4.5	30
	PrimaConnex™ WD 5.0, Genesis™ 5.5, 6.5	30
	Restore™ SD 3.3, 3.4	30
	Restore™ RD 3.75, 4.0	30
	Restore™ WD 5.0, 5.5, 6.0	30
<b>Nobel Biocare</b>	NobelActive™ NP 3.5	35
	NobelActive™ RP 4.3, 5.0	35
	NobelReplace™ NP 3.5	35
	NobelReplace™ RP 4.3	35
	NobelReplace™ WP 5.0	35
	NobelReplace™ 6.0	35
	Brånemark System™ NP 3.3	35
	Brånemark System™ RP 3.75, 4.0	35
	Brånemark System™ WP 5.0, 5.5	35
	SteriOss Replace™ 3.1, 3.5, 4.3, 5.0, 6.0	35
	SteriOss HL™ 3.25, 3.8, 4.5, 5.0, 6.0	35
<b>Straumann</b>	Bone Level 3.3 NC (Narrow CrossFit™)	35
	Bone Level 4.1, 4.8 RC (Regular CrossFit™)	35
	Standard/Standard Plus 4.8 RN (Regular Neck)	35
	Standard/Standard Plus 4.8 WN (Wide Neck)	35
	Standard Plus 3.3 NN (Narrow Neck)	35
<b>Sybron Implant Solutions</b>	Endopore 3.5, 4.1, 5.0 - external	30
	Entegra 3.5, 4.1, 5.0 - external	30
<b>Zimmer Dental</b>	Tapered Screw-Vent™ 3.5, 4.5, 5.7 / Tapered Screw-Vent 3.5, 4.5, 5.7 friction fit	30
	Screw Vent™ 3.3, 3.7, 4.5 / Screw-Vent 3.3, 3.7, 4.5 friction fit	30
	SwissPlus™ 4.8	35

**Note:** To remain aligned with the corresponding manufacturer's recommendations, torque values are subject to change without written notice.

All trademarks and company names are the property of their respective owners.

# We have you covered

Available for all major implant systems, ATLANTIS patient-specific prosthetic solutions are backed by comprehensive warranty terms\* for your added peace of mind.



All ATLANTIS cases are delivered with warranty tags as confirmation that an original ATLANTIS prosthetic solution was used for your case.

\*Subject to the terms and conditions of the ATLANTIS™ warranty.

## Medical and legal disclaimer

The information (including torque recommendations) is provided for information purposes only and does not constitute advice for medical or treatment purposes, nor a legal contract, between DENTSPLY Implants and any person or entity. Although all reasonable effort is made to present accurate information (including torque recommendations), DENTSPLY Implants makes no guarantees of any kind, including but not limited to that it may be outdated, invalid or subject to debate. The information is provided "as is" without warranty of any kind, either express or implied, including without limitation, the implied warranties of merchantability, fitness for use of a particular purpose, or non-infringement.

You agree that neither DENTSPLY Implants, nor any other company within the DENTSPLY international group, will be liable to you for any loss or damages, either actual or consequential, arising out of or relating to these terms, or to your (or any third party's) reliance upon information (including torque recommendations) obtained from DENTSPLY Implants. In particular, DENTSPLY Implants will have no liability for any consequential, indirect, punitive, special or incidental damages, whether foreseeable or unforeseeable, arising out of or relating to these terms, or to your (or any third party's) reliance upon information (including torque recommendations) obtained from DENTSPLY Implants, whether based in contract, tort, statutory or other law, except only in the case of death or personal injury where and only to the extent that applicable law requires such liability.